



# hekaDrill




## Ultimate Precision

*The hekaDrill has been developed with innovation and precision as a guiding principle.*

*Key capabilities include:*

- High speed motor control - up to 85,000 rpm. Allows high torque/low speed control. (min. 200 rpm)
- Advanced motor control ensures comfort for the surgeon during extended use with minimal heat gain.
- Adjustable bur exposure lengths.
- Wireless foot pedal operates using a secure encrypted network, ensuring greater user protection.
- Adjustable acceleration and braking functionality.



*A single rotary bur allows 5 possible exposure positions offering total flexibility and reduced inventory costs*

# ELECTRIC HANDPIECE AND ATTACHMENTS

HANDPIECE	Model	Order code	Speed rpm
-----------	-------	------------	-----------

HDS-121  
Cable length 4m



ATTACHMENTS	Model	Order code	Length - mm	Bur size - mm
-------------	-------	------------	-------------	---------------

## STRAIGHT

	S-C07	103360	70	2.38
	S-C08	103371	80	2.38
	S-C10	103383	100	2.38
	S-C15	103394	150	2.38
	S-D09	103409	90	3.18
	S-D14	103414	140	3.18

## ANGLED

	A-C07	103446	70	2.38
	A-C10	103451	100	2.38
	A-C15	103467	150	2.38
	A-D14	103472	140	3.18
	A-C07S	104070	70	2.28
A-C10S	104081	100	2.28	
A-C15S	104097	150	2.28	

## CRANIOTOME

	C-C09 Small	103508	90	2.38
	C-C10 Standard	103480	100	2.38
	C-D12 Large	103495	120	3.18

PERFORATOR ATTACHMENT	Model	Order code	Speed rpm	Chuck type
-----------------------	-------	------------	-----------	------------

P-P01 103519 1,000 Hudson



Used with HDS-121 at 72,000 rpm



Special order only

Locking mechanism spanner (for HDS-12x) Order code 103686 supplied with handpiece.





Description	Model	Order code
Console - less irrigation	PC-120	103305
Console - with irrigation	PC-121	103318



The advanced hekaDrill power console is available with or without irrigation and offers numerous features:

- Simultaneous dual handpiece operation.
- Adjustable acceleration and braking functionality.
- Integrated and customisable irrigation control with anti drip mode.
- 200+ user profiles can be stored and easily accessed.
- Rings illuminate to indicate port status for ease of use.

# IRRIGATION AND FOOT PEDALS

## STERILE IRRIGATION KIT



Sterile irrigation kit.  
Box of 5.  
Order code: S101461

## FOOT PEDALS



- Both wired and wireless foot pedal available.
- Wireless foot pedal operates using a secure 128 bit encrypted network for complete protection.
- A or B tool selection.
- Foot pedal underlighting to allow for easy detection.



4 metre cable for connecting and charging foot pedals via power console.

Order code: 103597

	Model	Order code
Wired	MFP-120	103556
Wireless	MFP-121	103562

## POWER SUPPLY

Mains power cable 3 metres long. Please specify from below.

UK	Order code 8290
Euro	Order code 8300
USA/Canada	Order code 8310
Australia/New Zealand	Order code 8320
Japan	Order code 8330
Swiss	Order code 14700



HANDPIECE	Model	Order code	Speed rpm
-----------	-------	------------	-----------



HDS-120	103339	100,000
---------	--------	---------

The HDS-120 handpiece accepts all straight and angled attachments and craniotome

PERFORATOR	Model	Order code	Speed rpm	Operating pressure
------------	-------	------------	-----------	--------------------



P-P02	103524	850	7 bar
-------	--------	-----	-------

Uses Hudson chuck type

FOOT CONTROL	Model	Order code	Foot pedal inlet connection	Foot pedal outlet connection
--------------	-------	------------	-----------------------------	------------------------------



FV-220	14900
--------	-------

Integral 3 metre hose with MA7 probe

Accepts 'wash assist' hose from sizes below



3 metre length - Order code 632333  
 4 metre length - Order code 632343  
 5 metre length - Order code 632353



Wash cap for use with HDS-120/ P-P02  
 Order code 9950



# STERILISATION TRAYS AND CONTAINER

## STERILISATION TRAYS AND CONTAINER



### ELECTRIC INSTRUMENT SYSTEM

Order code 103610  
Includes lid, container, electric tray and attachment tray.

Handpiece tray

Optional tray contents -

- 1 off HDS-121 handpiece
- 1 off P-P01 perforator
- 1 off electric cable
- 1 off spanner

Attachment tray

Optional tray contents -

- 5 off attachments
- 1 off spanner

### PNEUMATIC INSTRUMENT SYSTEM

Order code 103629  
Includes lid, container, pneumatic tray and attachment tray.

Handpiece tray

Optional tray contents -

- 1 off HDS-120 handpiece
- 1 off P-P02 perforator
- 1 off pneumatic hose
- 1 off spanner

Attachment tray

Optional tray contents -

- 5 off attachments
- 1 off spanner

### CONTAINER AND TRAYS ONLY

Wire container and lid

Order code 103944

Tray for electric system

Order code 103953

Tray for pneumatic system

Order code 103989

Tray for attachments

Order code 103966

UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)  
De Soutter Medical Limited  
Halton Brook Business Park, Weston Road  
Aston Clinton, Aylesbury, Bucks HP22 5WF  
Tel: +44 (0) 1296 634 000  
Fax: +44 (0) 1296 634 033  
Email: info@de-soutter.com

AUSTRALIA  
De Soutter Medical Australia Pty Ltd  
2/12-14 Apollo Drive, Hallam  
Victoria 3803  
Tel: +61 (0) 3 9702 4441  
Fax: +61 (0) 3 9702 4484  
Email: australia@de-soutter.com

AUSTRIA  
De Soutter Medical Austria  
Zweigniederlassung Österreich  
Dietrichsteingasse 10, A-3400 Klosterneuburg  
Tel: +43 (0) 676 96 71 770  
Fax: +43 (0) 2243 21 656  
Email: austria@de-soutter.com

BELGIUM  
De Soutter Medical Belgium  
Bessemersstraat 14  
3620 Lanaken  
Tel: +32 (0) 89/47 15 37  
Fax: +32 (0) 89/70 12 19  
Email: belgium@de-soutter.com

FRANCE  
De Soutter Medical France  
1 252 Avenue du Parc des Expositions  
33260 La Teste de Buch  
Tel: +33 (0) 5 56 54 89 36  
Fax: +33 (0) 9 70 61 37 60  
Email: france@de-soutter.com

GERMANY  
De Soutter Medical Deutschland  
Bahnhofstraße 4  
D-66625 Nohfelden  
Tel: +49 (0) 68 52-99 12 46  
Fax: +49 (0) 68 52-99 12 47  
Email: deutschland@de-soutter.com

ITALY  
De Soutter Medical Italia  
Località Fornace SNC  
27022 Casorate Primo - PV  
Tel: +39 (0) 2 9009 4098  
Fax: +39 (0) 2 9009 2673  
Email: italy@de-soutter.com

NETHERLANDS  
De Soutter Medical Nederland  
Gelderlandhaven 2X  
3433 PG Nieuwegein  
Tel: +31 (0) 85 0491480  
Fax: +31 (0) 85 0491489  
Email: nederland@de-soutter.com

UNITED STATES OF AMERICA  
De Soutter Medical USA, Inc.  
224 Rolling Hill Road, Suite 12A  
Mooresville, NC 28117  
Tel: +1 (704) 655 9040  
fax: +1 (704) 987 2035  
Email: usa@de-soutter.com

**deSoutter**  
MEDICAL

Distributed by De Soutter Medical, manufactured by Zethon